



Cardiovascular

One Scimed Place
Maple Grove, MN 55311-1566
763-494-1700 Tel
www.bostonscientific.com

TO: Medical Healthcare Professional
FROM: Boston Scientific Corporation
SUBJECT: Nickel Content - TAXUS[®] Liberté[®] Atom[™] Stent System
DATE: 21-May-2009

Per your request for additional information, the TAXUS[®] Liberté[®] Atom[™] Stent is made from 316L surgical grade stainless steel tubing, and consists of the following elements:

| Element | Elemental Composition By Weight ¹ |
|------------|--|
| Iron | 65% |
| Nickel | 14% |
| Chromium | 18% |
| Molybdenum | 3% |

Section 4 (Warnings) of the Directions for Use

“Patients with known hypersensitivity to 316L stainless steel may suffer an allergic reaction to this implant.”

To receive a complete copy of the DFU, please visit www.bostonscientific.com.

¹ (Source: American Society for Testing of Materials (ASTM) Standards: ASTM F138-03 (316L Stainless Steel))